


Workflow Improvements from SARSTEDT to Reduce ToTAT and Improve Patient Care

Reduce recollections with the S-Monovette® blood collection system

Better for the lab, better for the patient

- Provides the option for gentle syringe-style collection directly into treated, analyzer-ready tubes
- Eliminates sample transfers from untreated syringes
- Reduces hemolysis and clotting rates
- Ultra-low volumes for better blood management



Contact us for a lunch-and-learn or automation needs analysis



Streamline Transport with the Tempus600® Specimen Transportation System

It's like having a lab right next to the patient

- Dedicated system for specimen tubes
- Tubes placed directly into sending stations - no packing or batching required
- Predictable delivery to the laboratory within seconds
- Evens out peaks for continuous flow of samples to support FIFO

Integrate Pre- and Post-Analytical Automation

Let us create an independent or cooperative solution for you

- Functional options include: bulk or rack sample loading; accessioning; sorting to racks or bins; decapping; aliquoting; transfer to laboratory line; recapping for archival; recapping for send-out testing
- Combine with Tempus600® systems and samples are racked and ready for analysis within minutes after collection



For more information about pre- and post-analytical solutions from SARSTEDT, or for a free needs assessment, see the link below.

<https://www.sarstedtsamples/Workflow.html>